

# *Aids for People with Visual or Physical Disabilities*

<b>Product Name (Manufacturer)</b>	<b>Description/Use</b>
<b>Syringe and Pump Magnifiers</b>	
<b>BD Magni-Guide</b> (BD)	Syringe is slipped into curved channel of Magni-Guide. At opposite end, insulin vial is snapped into collar. Guide magnifies entire scale 1.7 times. Uses BD 1-cc, ½-cc, or ⅜-cc syringe. Available from BD, local pharmacies, Medicool, and Gemco Medical.
<b>Pump Monocle</b> (Magnifics)	Display magnifier for insulin pumps. Two styles. Soft Monocle wraps around the device with 2 to 3 times magnifying lens; or flip version.
<b>Syringe Magnifier</b> (Apothecary Products)	Clips firmly to syringe barrel. Easy calibration readings of ½-cc and 1-cc syringes. Uses any standard-size syringe.
<b>Tru-Hand</b> (Whittier Medical)	Insulin bottle holder/syringe guide/magnifier. Locks in syringes with direct guide into vial stopper. Magnifies entire scale. Insulin bottles can be changed without removing syringe from device.
<b>Needle Guides, Vial Stabilizers, and Vial Syringe Guides</b>	
<b>Inject Assist</b> (Apothecary Products)	Holds syringe and insulin bottle steady as syringe is filled. Preset dosage guide for accuracy.
<b>Injection Safety Guard</b> (Apothecary Products)	Protects hands from accidental injection. Fits over the top of insulin vial, allowing user to place hands underneath safety guard and begin filling the syringe.
<b>NeedleAid</b> (NeedleAid)	Stabilizing guide for injection using an insulin pen or syringe. The broad base indicates where the injection will be made. Ensures injection at the proper angle and to proper depth. Shields the needle for increased safety. Suitable for use with one hand. Insulin is injected manually. Needle is hidden from view to reduce needle fear and pain perception.
<b>Vial Syringe Guide</b> (Apothecary Products)	Eliminates blunted or broken needles by safely connecting syringe to insulin bottle.