

# Insulin Pens

Product Name (Manufacturer)	Comments
<b>Autopen</b> (Owen Mumford)	Reusable pen for use with Eli Lilly insulin cartridges. Available in four models. AN3000 (1.5-ml insulin cartridge) delivers 2 to 32 units in 2-unit increments. AN3100 (1.5-ml cartridge) delivers 1 to 16 units in 1-unit increments. AN3800 (3-ml cartridge) delivers 2 to 42 units in 2-unit increments. AN3810 (3-ml cartridge) delivers 1 to 21 units in 1-unit increments. Side button releases dose. Dosage indicator returns to start line after injection.
<b>Humalog KwikPen</b> (Eli Lilly)	Prefilled disposable pen for use with either Humalog, Humalog Mix 75/25, or Humalog Mix 50/50. Holds 3 ml of insulin and delivers 1 to 60 units in 1-unit increments. Fix over-dialing by dialing knob backward. Refrigerate unused pens. In-use pen should be stored at room temperature.
<b>Original Humalog Pen, Humalog Mix 75/25 Pen, Humalog Mix 50/50 Pen, Humulin N Pen, and Humulin 70/30 Pen</b> (Eli Lilly)	Disposable pens prefilled with Humalog, Humalog Mix 75/25, or Humalog Mix 50/50. Holds 3 ml of insulin and delivers 1 to 60 units in 1-unit increments. Audible clicks while dialing. Fix over-dialing by dialing knob backward. Refrigerate unused pens. In-use pen should be stored at room temperature.
<b>HumaPen Luxura HD</b> (Eli Lilly)	Reusable pen for use with Humalog insulin. For children with diabetes. Holds a 3-ml cartridge and delivers 1 to 30 units in ½-unit increments after first unit. Fix over-dialing by dialing knob backward. Do not store pen in refrigerator.
<b>HumaPen Memoir</b> (Eli Lilly)	Reusable pen for use with Humalog insulin. With memory feature that digitally displays last 16 doses. Holds a 3-ml cartridge and delivers 1 to 60 units in 1-unit increments. Fix over-dialing by dialing knob backward. Do not store pen in refrigerator. Store unopened insulin cartridges in refrigerator.
<b>Lantus SoloStar</b> (Sanofi-Aventis)	Disposable pen prefilled with Lantus. Holds 3 ml of insulin and delivers 1 to 80 units in 1-unit increments. Dosage indicator returns to 0 after injection. Refrigerate unused pens. In-use pen should be stored at room temperature.
<b>Levemir FlexPen</b> (Novo Nordisk)	Disposable pen prefilled with Levemir. Holds 3 ml of insulin and delivers 1 to 60 units in 1-unit increments. Audible clicks while dialing. Dosage indicator returns to 0 after injection. Fix over-dialing by dialing knob backward. Refrigerate unused pens. In-use pen should be stored at room temperature.
<b>NovoLog FlexPen</b> (Novo Nordisk)	Disposable pen prefilled with NovoLog. Holds 3 ml of insulin and delivers 1 to 60 units in 1-unit increments. Audible clicks while dialing. Dosage indicator returns to 0 after injection. Fix over-dosing by dialing knob backward. Refrigerate unused pens. In-use pen should be stored at room temperature.
<b>NovoLog Mix 70/30 FlexPen</b> (Novo Nordisk)	Disposable pen prefilled with NovoLog Mix 70/30. Holds 3 ml of insulin and delivers 1 to 60 units in 1-unit increments. Audible clicks while dialing. Dosage indicator returns to 0 after injection. Fix over-dialing by dialing knob backward. Refrigerate unused pens. In-use pen should be stored at room temperature.
<b>NovoPen Junior</b> (Novo Nordisk)	Pen holds a 3-ml Novo Nordisk PenFill cartridge (use with NovoLog and Novolin) and delivers 1 to 35 units in ½-unit increments. Dosage indicator returns to 0 after injection. Refrigerate unused cartridges. In-use pen should be stored at room temperature.
<b>NovoPen 3</b> (Novo Nordisk)	Disposable pen for use with NovoLog and Novolin insulins. Holds 3 ml of insulin and delivers 2 to 70 units in 1-unit increments. Dosage indicator returns to 0 after injection. Do not refrigerate.
<b>OptiClik</b> (Sanofi-Aventis)	Reusable pen for use with Lantus or Apidra. Holds 3-ml cartridge and delivers 1 to 80 units in 1-unit increments. Cartridge clicks into place to ensure proper setup. Dose displayed digitally. Dosage indicator returns to 0 after injection. Do not refrigerate pen. Store unopened insulin cartridges in refrigerator.